

SCHOTT TopLy[®] Vials



General Product Information

SCHOTT TopLy[®] vials consist of SCHOTT Type I glass combined with a transparent and nonporous Si-O-C-H layer. They comply with all current standards, such as Ph.Eur., USP and JP.

Due to the coating along with an optimized container geometry, this vial is especially recommended for lyophilization. SCHOTT TopLy[®] offers the advantage of less rejects and minimized residual volume after reconstitution. Furthermore, the removal of the drug is improved.

Physical & Chemical Properties

The ultra thin layer of SCHOTT TopLy[®] vials is characterized by the following properties:

Physical Data
Layer thickness of ~ 40 nm
Stable against mechanical load
Stable washing process
Stable sterilization: <ul style="list-style-type: none"> · Autoclaving (121 °C) · Depyrogenation (dry heat treatment at 250 °C – 350 °C)

Chemical Data
Chemical layer properties: Si-O-C-H
Long-term stable layer system during storage proven by accelerated aging at 40 °C and room temperature after 3 years
Bond covalently to the material and chemically uniform
Dense coating (i.e. non porous)
Contact angle for water: ~100° (hydrophobic surface)

Verifications


Extractable volume

Method:

- Take out 4 ml solution from 10 ml vials containing 12.5 mg/ml IgG1 and 12.5 mg/ml mannitol dissolved in a histidine/glycine mixture
- Buffered to pH 6.0 with a needle through the stopper
- Extractable volume determined gravimetrically (n=3)

Result:

- ~4% more extractable volume from hydrophobic SCHOTT TopLy[®] compared to standard Type I glass
- 4% equals 160 µl and contains 2 mg of protein




Contact angles measured over 36 months at room temperature

Method:

- 1400 SCHOTT TopLy[®] vials, 2 ml
- WFI, phosphate and citrate buffer (pH 5, 7, 8)
- Sterilization at 121 °C, 30 min
- Samples stored at 40 °C


Result:

- Contact angle > 90°



Product Information

Thanks to our patented coating technology, vials are endowed with a stable and covalently bond hydrophobic interior. This coating promotes efficient liquid removal and enhanced drug stability.



Maximum inspection – validated process

Stage 1

Two 100% on-site inspections on each reactor (temperature, optical plasma emission)

Stage 2

Control of process parameters (on-line, including gas flow, vacuum, microwaves)

Stage 3

Automatic System Monitoring of long-term stability (maintenance, calibration of the actuators and sensor, data acquisition and long term storage)

Value-adding Product Benefits and Services

Application ranges

Improved lyo cake

Less disruption and production losses due to reduced climbing at vial edges

Optimized yield

Reduced adsorption at the inner glass surface minimizes the residual volume after reconstitution

High process stability

Diminished glass breakage as stress risk is low due to special geometry design

Process time reduction

Optimized bottom geometry allows for better heat transfer during the freeze drying process

Dimensional changes not required

Due to the very thin coating layer (~40 nm), no dimensional changes are needed

Packaging

- SCHOTT TopLyo® vials are delivered in special trays with optional separators to avoid glass to glass contact
- A standard Euro Pallet (1200 x 800 mm) contains 15–27 layers of 9 trays each

Capacity	2 R	6 R	10 R
Pieces/tray	344	186	154