



## Stock Rubber Stoppers



We hold a large range of rubber stoppers available for immediate dispatch from our warehouse

• Injection stoppers • Freeze Dry stoppers • Diagnostic stoppers •

INJECTION STOPPERS / INFUSION STOPPERS												
Adelphi Part No.	Description	West SAP Code	West Art.	Formulation	Final treatment	Lamination	Silicone	Packaging	Pack Qty	Carton Qty	MOQ	
INJ13TB3WRSNP1	13mm Injection Stopper FluroTec®	7002-5046	1358	4023/50 Grey	NovaPure® RS	Daikyo fluro resin D	B2-40	5 x STERIlizableBAG™	1,000	5,000	5,000	
INJ20TB3WRSNP	20mm Injection Stopper FluroTec®	7002-5047	1343	4023/50 Grey	NovaPure® RS	Daikyo fluro resin D	B2-40	2 x STERIlizableBAG™	3,000	6,000	6,000	
INJ13TB3WSRS	13mm Injection Stopper FluroTec®	1012-2737	1358	4023/50 Grey	Westar Select® RS	Daikyo fluro resin D	B2-40	STERIlizableBAG™	2,000	10,000	2,000	
INJ20TB3WSRS	20mm Injection Stopper FluroTec®	1014-5401	1343	4023/50 Grey	Westar Select® RS	Daikyo fluro resin D	B2-40	STERIlizableBAG™	2,000	6,000	2,000	
INJ13TBWRS	13mm Injection Stopper Flurotec®	7001-5950	1258	4110/40 Grey	Westar® RS	Flurotec® +	B2-40	STERIlizableBAG™	5,000	20,000	Please Contact Us	
INJ20TBWRS	20mm Injection Stopper Flurotec®	7001-4383	1243	4110/40 Grey	Westar® RS	Flurotec® +	B2-40	STERIlizableBAG™	1,000	5,000	1,000	
1970-0004	13mm Injection Stopper FluroTec®	1970-0004	S2-F451	4432/50 Grey	Westar® RS	Daikyo fluro resin D	B2-40	STERIlizableBAG™	1,000	5,000	1,000	
INJW13RtS	13mm Injection Stopper	7001-5250	1104	4023/50 Grey	RTS™	none	A	STERIlizableBAG™	4,000	16,000	4,000	
INJW20RtS	20mm Injection Stopper	7001-7999	1071	4023/50 Grey	RTS™	none	A	STERIlizableBAG™	1,500	4,500	1,500	
INJW20S	20mm Injection Stopper Silicone	7000-4778	1185	8763/40	Standard	none	none	Foilbag	1,500	4,500	1,500	
INJ20WC	20mm Injection Stopper White	7000-3332	1213	4601/45 White	Standard	none	none	Polybag	7,000	7,000	500	
INJW32G3RTS	32mm Infusion Stopper	7002-1427	4440	4023/50 Grey	RTS™	none	A	STERIlizableBAG™	500	1,500	500	
FREEZE DRY STOPPERS												
FD13TT3WRS	13mm Freeze Dry Stopper FluroTec <sup>4</sup>	7001-9812	1356	4023/50 Grey	Westar® RS	Daikyo fluro resin D	none	STERIlizableBAG™	2,000	10,000	Please Contact Us	
FD13TTB3WRS	13mm Freeze Dry Stopper FluroTec <sup>4</sup>	7002-5173	1356	4023/50 Grey	Westar® RS	Daikyo fluro resin D	B2-TR	STERIlizableBAG™	2,000	10,000	2,000	
FD20TT3WRS	20mm Freeze Dry Stopper FluroTec <sup>4</sup>	7001-9820	1346	4023/50 Grey	Westar® RS	Daikyo fluro resin D	none	STERIlizableBAG™	1,000	4,000	Please Contact Us	
FD20TTB3WRS	20mm Freeze Dry Stopper FluroTec <sup>4</sup>	7002-5169	1346	4023/50 Grey	Westar® RS	Daikyo fluro resin D	B2-TR	STERIlizableBAG™	2,000	4,000	2,000	
FD13TBrtS	13mm Freeze Dry Stopper Flurotec®	1956-0585	V2-F8	D713 Grey	RTS™	Daikyo fluro resin D	B2-40	STERIlizableBAG™	5,000	10,000	5,000	
FD20TBrtS	20mm Freeze Dry Stopper Flurotec®	1956-0589	V10-F22	D713 Grey	RTS™	Daikyo fluro resin D	B2-40	STERIlizableBAG™	1,750	3,500	1,750	
FDW13WRS	13mm Freeze Dry Stopper	7002-1585	1079	4023/50 Grey	Westar® RS	none	4.5	STERIlizableBAG™	3,000	12,000	3,000	
FDW20RtS	20mm Freeze Dry Stopper	7000-5780	1319	4023/50 Grey	RTS™	none	A	STERIlizableBAG™	2,000	4,000	2,000	
DIAGNOSTIC STOPPERS												
FDIA14WGBRTS	14mm Freeze Dry Stopper	7001-1322	1067	W1816 Grey	Standard	none	B	STERIlizableBAG™	3,000	9,000	3,000	
FDIA18W/B	18mm Freeze Dry Stopper	7000-5178	1008	PH701/40 Red	Standard	none	B	Polybag	7,000	7,000	500	
FDIA18WG/B	18mm Freeze Dry Stopper	7001-1079	1008	W1816 Grey	Standard	none	B	Polybag	7,000	7,000	500	
FDIA22	22mm Freeze Dry Stopper	7000-1692	1010	PH701/40 Red	Standard	none	A	Polybag	3,000	3,000	500	
SDIA18RtS	18mm Plug Stopper	7000-5215	1007	PH701/40 Red	RTS™	none	A	STERIlizableBAG™	4,000	8,000	4,000	

### FluroTec® (FluroTec® + and Daikyo fluro resin D)

An ETFE Film applied to the stoppers to significantly reduce the risk of interaction of the drug product with the closure

### Ready to Sterilize (RTS™) and Westar® RS (WRS)

Final rinse in WFI and packed in class 100 (ISO 5 clean room into a STERIlizableBAG

### B2-40

Cross linked silicone coating applied to the top surface of the rubber stopper

Adelphi Healthcare Packaging  
Olympus House, Mill Green Road  
Haywards Heath, West Sussex  
RH16 1XQ, UK  
Tel: +44 (0)1444 472300  
Email: sales@adelphi-hp.com  
Web: www.adelphi-hp.com