

True to Our Core™



## Westar® RU Syringe Plungers



### A Ready-to-Use, Steam-Sterilized Solution for Prefillable Glass Syringes

Westar® Ready-to-Use (RU) plungers can mitigate risk for chemical and physical deterioration of the component through an optimized autoclave steam sterilization cycle that, when compared to high-dose gamma irradiation, may:

- Minimize impact on the elastomer physical characteristics, and on chemical oxidation and degradation during the plungers' shelf-life
- Decrease the number and lower the level of extractables, with less potential for interaction between the plunger and the drug product
- Improve prefilled syringe functionality due to lower and more consistent break loose forces

Plungers in primary configurations are available for rapid order fulfillment. West and its partner, Daikyo Seiko, Ltd., offer primary configurations in market-proven elastomer formulations. Westar RU steam sterilized plungers are available with or without FluroTec® film, fit standard prefillable syringes and comply with USP, JP and Ph. Eur. Pharmacopeia requirements.

The Westar RU steam sterilization process is fully validated and backed by the technical and regulatory expertise and support that West is known for.

#### FOR MORE INFORMATION

**West Pharmaceutical Services, Inc.**  
North America: +1 800 345 9800  
Europe: +49 2403 79 60  
Asia Pacific: +65 6860 5879  
[www.WestPFSolutions.com](http://www.WestPFSolutions.com)  
[www.westpharma.com](http://www.westpharma.com)