



## adaptiQ® Sterile Injection Vials

SCHOTT ISO standard sterile and depyrogenated, ready-to-use vials, supplied in a choice of three packaging formats

### SCHOTT adaptiQ® Vials in Cup Nest

SCHOTT Injection Vials supplied sterile nestled in cups within a plastic preform, inside an industry standard tub. Cup Nest packaging offers simple de-nesting and re-nesting and is also ideal for automated processing.



### SCHOTT adaptiQ® Vials in Tray

Sterile vials in a tray ideal for large production runs using traditional infeed equipment, or for situations where full removal of the vials is required for manual filling.



### SCHOTT adaptiQ® Vials in Clip Nest

Featuring a nest with clips holding each sterile vial securely by the body, clip nest packaging is compatible with fully-nested aseptic, automated processing, and is also suitable for freeze drying in the nest.



Cup Nest NBB (No Blowback)				
Size	Adelphi Part No.	Quantity per Tub	MOQ	Quality
2R	VC2RNBBBQC	100	15 Tubs	FHU
6R	VC6RNBBBQC	48	15 Tubs	FHU
10R	VC10RNBBBQC	48	15 Tubs	FHU
20R	VC20RNBBBQC	24	15 Tubs	FHU
50R	VC50RNBBBQC	16	15 Tubs	FHU

Tray NBB (No Blowback)				
Size	Adelphi Part No.	Quantity per Tray	MOQ	Quality
2R	VC2RNBBBQ	160	10 Trays	FHU
6R	VC6RNBBBQ	88	10 Trays	FHU
10R	VC10RNBBBQ	80	8 Trays	FHU
20R	VC20RNBBBQ	48	8 Trays	FHU
50R	VC50RNBBBQ	28	6 Trays	FHU

Clip Nest NBB (No Blowback)				
Size	Adelphi Part No.	Quantity per Tub	MOQ	Quality
2R	VC2RNBBBQ	100	15 Tubs	FHU
6R	VC6RNBBBQ	48	15 Tubs	FHU
10R	VC10RNBBBQ	48	15 Tubs	FHU

Clip Nest EBB (European Blowback)				
Size	Adelphi Part No.	Quantity per Tub	MOQ	Quality
2R	VC2REBBBQ	100	15 Tubs	FHU
10R	VC10REBBBQ	48	15 Tubs	FHU

Unsure which secondary packaging option is best suited to your process?  
Contact our friendly sales team and we will be happy to assist you:

[sales@adelphi-hp.com](mailto:sales@adelphi-hp.com)