

## MixJect

### Integrated transfer device for drug reconstitution and injection

The MixJect® system enables the safe, rapid and easy preparation of lyophilized drugs. The system is a single unit for reconstituting a powder drug with a diluent prefilled syringe. Upon reconstitution, the drug is available for immediate injection with a pre-attached needle.

#### MixJect System Attributes

- Preamsembled with a client-specified standard or safety needle
- Fits standard vials, syringes, needles and Luer lock connections
- Compatible with popular auto-injector devices
- Designed for minimal residual volume and overflow



#### MixJect is ideal for

- Lyophilized drugs supplied with prefilled diluent syringes
- Depot injections
- Self administration of drugs for chronic illness
- The MixJect system is 510(k) cleared by the United States Food and Drug Administration and carries the CE mark.

#### Innovations in Safety and Administration Systems

West's administration systems enhance the safe delivery of drugs through advanced reconstitution, mixing and transfer technologies. The systems are easy to use, help promote safety and compliance and use your existing formulation, manufacturing process and primary packaging. West has extensive experience helping our customers around the world navigate local regulatory requirements to facilitate registration.

To find out how West's innovative systems can add value to your drug product or application, please contact your local West account representative.

For more information visit [www.westpharma.com](http://www.westpharma.com)

**North America** 1-800-345-9800 option 8 | **South America** +55 11 4055 6061 | **Europe** +49 (0) 2403-7960 | **Asia Pacific** +65 6862 3400

**West Pharmaceutical Services, Inc.** 530 Herman O. West Drive | Exton, PA 19341